




MEDRUX®

PRODUCT GUIDE

Yummylooks Sdn Bhd



MEDRUX®

Company Name:

Yummylooks Sdn Bhd

Address

Unit 27, 3A-01, Glo Damansara
699 Jalan Damansara
TTDI, Kuala Lumpur

Phone

+6 03 – 29359191

Whatsapp

+6 017 279 5502

Email / Web

client@medrux.com

www.medrux.com



Product Name:
ATOMX

Product Type:
Sterile Latex Surgical Gloves

1. DESCRIPTION AND MATERIALS

Product : Gamma Ray / EtO Sterilized Powder Free Latex Surgical Gloves, Hand specific; curved fingers; natural; textured surface; USP grade absorbable corn starch

Materials : 100% Rubber Latex

Cuff : Beaded

Sizes : Extra-Small; Small; Medium; Large and Extra-Large

Colour : White

Compliance: CE, ASTM D3577, EN455 , ISO13485



Product

ATOMX - Sterile Latex Surgical Gloves




2. PERFORMANCE REQUIREMENTS

Gloves are sampled at random and inspected in accordance with ISO2859-1:1999 based on batch size. The inspection levels and Acceptance Quality Levels (AQL) shall conform to specification as shown below or as agreed between the purchaser and the supplier.

Samples are inspected/tested as per ASTM D3577 and EN 455 Standards

Storage	The gloves shall maintain their properties when stored in a dry condition at temperature not higher than 30° C.
Moisture Content	Below 0.8% per glove.
Shelf-life	5 years from the date of manufacturing.

3. DIMENSIONS

Dimension	Standards		
	MEDRUX	ASTM D3577	EN455
Length(mm)	Min 280	Min 245 (size 5.5) Min 265 (size 6-9)	Min. 250 (size 5-5.5) Min. 260 size 6-6.5) Min. 270 (size 7-8) Min. 280 (size 8.5-9.5)
Width(mm) • 5.5 • 6.0 • 6.5 • 7.0 • 7.5 • 8.0 • 8.5 • 9.0	71 ± 5 77 ± 5 83 ± 5 89 ± 5 95 ± 5 102 ± 5 108 ± 5 115 ± 5	70 ± 6 76 ± 6 83 ± 6 89 ± 6 95 ± 6 102 ± 6 108 ± 6 114 ± 6	72 ± 4 77 ± 5 83 ± 5 89 ± 5 95 ± 5 102 ± 6 108 ± 6 114 ± 6
Thickness - single wall (mm)  Fingers  Palm  Cuff	Min 0.10 Min 0.10 Min 0.10	Min 0.10 Min 0.10 Min 0.10	N/A

Product Type:

ATOMX - Sterile Latex Surgical Gloves

4. PHYSICAL PROPERTIES

Size: All sizes

Property	ASTM D3577	EN455
Elongation At Break% ● Before Aging ● After Aging	min 750 min 560	N/A
Tensile Strength(MPa) ● Before Aging ● After Aging	min 24 min 18	N/A
Force at Break (N) ● Before Aging ● After Aging	N/A	Median 9 Median 9

5. FREEDOM FROM HOLES

Gloves shall meet requirements as stipulated in ASTM D3577 and EN 455.

6. POWDER RESIDUE

Powder Residue is tested as per ASTM D 6124-06. HLRI Powder residue shall be > 2mg/glove.

7. FUNCTIONAL BENEFITS

- *Extra strength provides additional protection from surgical debris.*
- *Fully anatomical design to reduce hand fatigue.*
- *Softness provides superior comfort and natural fit.*
- *Micro-roughened surface provides excellent wet and dry grip.*
- *Easy donning and helps prevent rolling back.*
- *High strength and elasticity.*

Product Type:

ATOMX - Sterile Latex Surgical Gloves

8. PACKAGING

Pouch: 1 pair


Inner Box: 50 pairs

Outer Carton: 200 pairs (4 inner boxes)

9. MARKING

The dispenser and outer cartons shall have the following information/prints:

- *Product description and Brand including size and net quality of gloves.*
- *Name and country of origin of manufacturer or Distributor.*
- *Single Use only or the symbol for symbol for single use.*
- *Other instructions on storage.*
- *Manufacturer's or Distributor's identification i.e. Lot Number; Bar Codes; Reorder Number.*



Product Photos



CE COMPLIANCE

Certificate of Compliance

This certificate of compliance is here by issued to the below applicant

Applicant — YUMMYLOOKS SDN BHD

Address — 27, 3A-01, Glo Damansara, TTDI, 6000, Kuala Lumpur Malaysia

Manufacturer — Medrux

Product — Latex and Nitrile Gloves

Trade Name YUMMYLOOKS SDN BHD

Applicable Directive(s) — Medical Devices Directive -93/42/EEC

Test standard — EN 455 Part 1:2000 / EN 455 Part 2:2015/ EN 455 Part 3:2015

Test Report no — CRSSA/200845255-CA45230

Test Laboratory name — SGS (MALAYSIA) SDN BHD

This Certificate of Compliance is issued under following conditions:

1 - The product is declared by the applicant as complying with the applicable directives.

2 - EA has performed a review of the technical information of the product.

3 - The applicant is required to inform EA of any changes to its product design and technical documentation.

4 - The CE mark may be used by the applicant, for this product, under the sole responsibility of the applicant*

Certificate Number: 21507

Initial Registration Date : 25/01/2021

Date of Certificate Issue : 28/01/2021

1st Surveillance Audit : 27/01/2022

2nd Surveillance Audit : 26/01/2023

Certificate Expiry* : 25/01/2024



Harshavardhan Reddy
(Director)

EURO ASSESSMENTS & CERTIFICATION LTD
Kemp House, 160 City Road, London, EC1V 2NX, United Kingdom

*Validity of the certificate is subjected to successful completion of surveillance audit on or before of due date (in case surveillance audit is not allowed to be conducted, this certificate shall be suspended/withdrawn).

Certificate Verification: Please verify certificate at <https://www.euroassessments.co.uk>. This certificate is the property of Euro Assessments & Certification Ltd and shall be returned immediately when demanded.

*The product liability rests with the manufacturer in accordance with their committed specifications and contractual obligations, including applicable standards. The CE mark shown above is for reference only and does not indicate accreditation. Euro Assessments does not accept any liability arising out of this certificate.





CERTIFICATE

This is to Certify that the Management System of

YUMMYLOOKS SDN BHD

3A-01A, Glo Damansara TTDI, Kuala Lumpur, Malaysia

has been found to conform to the Medical Devices – Quality Management System standard:

ISO 13485:2016

This certificate is valid for the following scope of operations:

“Latex and Nitrile Examination Gloves (Powdered and Powder-free)”

Certificate No.: MY91540H

Date of initial registration	Date of this Certificate	Surv. audit on or before/ Certificate expiry	Recertification Due
03 February 2021	03 February 2021	02 February 2022	02 February 2024

This Certificate remains valid subject to satisfactory surveillance audits.




Director

For verification and updated information concerning the present certificate visit to www.lmscert.com. This Certificate is the property of LMS Certification Limited and shall be returned immediately when demanded.



LMS Certification Limited
Labrynth Business Centre, 43 Middle Hill Gate, Stockport, Great Manchester, England-SK1 3DG
Phone :+44 208 935 5094
Company No.: 11029176
Visit :- www.lmscert.com
E-mail :- info@lmscert.com

LMSFM/01/MDQ/REV06

MEDRUX®

8



SFDA APPROVED FOR LATEX POWDER FREE

الهيئة العامة للغذاء والدواء

Saudi Food & Drug Authority

SFDA

About SFDA

Information Lists

Areas

Consumer Corner

Media Centre

Eservices

Search

Product Name

Manufacturer Name

low Risk Device Number

Product Number

Account Name	No value
Low risk device number	LR-20213-C-3362--SFDA-1037
model Number	MX9965 (Small, Medium, Large, X-Large)
Manufacturer Name	Yummylooks Sdn Bhd
Product name	Medrux Latex Gloves Powder Free
product Number	LR-20213-C-3362--SFDA-1037
Expired Status	Not active
Suspension Cancelation Status	Active
Manufacturer Country	Malaysia

مصنع

مها طالب عبدالله السعدون للصناعة

MTS II

LR-20193-C-0212-SFDA-0003


مصنع

مها طالب عبدالله السعدون للصناعة

MTS III

LR-20193-C-0212-SFDA-0004

[Check Here](#)



Business Registration & Rubber License



SURUHJAJAYA SYARIAH MALAYSIA
COMPANIES COMMISSION OF MALAYSIA
(Agensi di bawah KPDNKK)



COMPANIES ACT 2016
(ACT 777)

[Section 17]

**CERTIFICATE OF INCORPORATION
OF PRIVATE COMPANY**

This is to certify that

YUMMYLOOKS SDN. BHD.
(1245512-W)

is, on and from the 05th day of September 2017, incorporated under the Companies Act 2016, and that the company is a company limited by shares and that the company is a private company.

Dated at KUALA LUMPUR this 05th day of September 2017.



DATO' ZAHRAH ABD WAHAB FENNER
REGISTRAR OF COMPANIES
MALAYSIA



CHE ADERI BIN CHE AHMAD
Licensed Company Secretary
LS 0008651

CERTIFIED TRUE COPY




Info

DocID : 335562 Date : Wed Sep 06 10:11:36 +08 2017
Printing Date : 06-09-2017

This certificate is generated from SSM e-Info Services on at 06-09-2017 10:11:36
MENAPU SUBSISTENSI, NO. 7 JALAN SENTRAL 5, KUALA LUMPUR SENTRAL, 50470 KUALA LUMPUR.
Tel: 63-2289 4400 Fax: 63-2289 5411

SALINAN PEMOHON



LEMBAGA GETAH MALAYSIA
MALAYSIAN RUBBER BOARD

LESEN
[Peraturan 5(1)]

NO. SIRI : 036553
LGM/BEPL/W/01/3/018/2020

**AKTA LEMBAGA GETAH MALAYSIA (PERBADANAN) 1996
PERATURAN-PERATURAN LEMBAGA GETAH MALAYSIA
(PELESENAN DAN PERMIT) 2014**

NO. LESEN : W/01/16671

Mengikut Peraturan-Peraturan Lembaga Getah Malaysia (Pelesenan dan Permit) 2014 dan tertakluk kepada syarat-syarat yang dikenakan pada lesen ini, Lembaga mengeluarkan lesen yang berikut:

Jenis Lesen	Aktiviti Yang Dibenarkan	Aktiviti Tambahan
Mengeksport Sarung Tangan Getah	Mengeksport Sarung Tangan Getah	


menjalankan aktiviti atas nama

YUMMYLOOKS SDN BHD (1245512 - W)

bagi tempoh 24/08/2020 hingga 23/08/2021 di premis yang berikut :

YUMMYLOOKS SDN BHD
3A-01, GLO DAMANSARA, 699, JALAN DAMAI
TAMAN TUN DR ISMAIL
60000 KUALA LUMPUR
WILAYAH PERSEKUTUAN

Catatan : Syarat-syarat yang dikenakan pada lesen ini adalah seperti berikut :
i, iii & xxvii



Amir Hassan
Pegawai Yang Diberi Kuasa

*Sila lihat sebelah



OTHER DETAILS

LOCATION

Malaysia

PRODUCTION DURATION

21- 30 DAYS

Quality Inspection

At Warehouse

PAYMENT

50% down payment & 50% after
inspection at warehouse

CERTIFICATIONS

CE / LGM / ISO13485/ ISO9001 / DOC

BRAND

MEDRUX



1. MEDRUX® issues an offer to buyer
2. Purchase Order (PO) from buyer to be presented to MEDRUX.
3. After buyer acceptance, MEDRUX issues sales and purchase agreement (SPA).
4. Buyer signs the (SPA)
5. MEDRUX issues Invoice.
6. Payment terms 50% TT upfront and 50% after quality inspection at the warehouse.
8. Shipping to be arranged by MEDRUX or buyer.
9. Shipment to port

Company Name:

Yummylooks Sdn Bhd

Address

Unit 27, 3A-01, Glo Damansara

699 Jalan Damansara

TTDI, Kuala Lumpur

Phone

+6 03 – 29359191

Whatsapp

+6 017 279 5502

Email / Web

client@medrux.com

www.medrux.com



MEDRUX®